# User Manual of EpiProfile

#### Garcia Lab at UPenn

(bgarci@pennmedicine.upenn.edu, zuoyuan@pennmedicine.upenn.edu) Ref: Molecular & Cellular Proteomics 2015, 14 (6): 1696-707

#### **Contents**

- A. Requirements
- **B.** Instructions
- C. Steps
- D. Trouble shooting

# A. Requirements:

- 1. Windows system (Windows 7 or later version, any Intel or AMD x86-64 processor, RAM with 2 GB or more).
- 2. Xcalibur and Matlab should be installed firstly.

#### **B. Instructions:**

- 1. modify the input parameters (open the folder 'EpiProfile', open the file 'paras.txt', put your data path after 'raw\_path', set other parameters following the instructions).
- 2. start Matlab (in the folder 'EpiProfile', doubly click the file 'EpiProfile.m').
- 3. run EpiProfile (in the Matlab Command Window input "EpiProfile" and press "Enter").

The results are under the data path (histone\_layouts and histone\_ratios.xls, or histone\_ratios\_SILAC.xls, or histone\_ratios\_C13.xls, or histone\_ratios\_N15.xls, or histone\_ratios\_13CD3.xls).

### C. Steps:

Step 1 (modify the input parameters)

### [EpiProfile]

% the datapath of raw files raw\_path=C:\F\Exp76\1

#### % 1: Human, 2: Mouse

norganism=1


% 1: histone\_normal, 2: histone\_SILAC, 3: histone\_C13, 4: histone\_N15, 5: histone\_13CD3

nsource=1


% if histone\_N15, 0: N14 light Mods, 1: N15 light Mods, 2: N14 heavy Mods, 3: N15 heavy Mods, 4: 0+1, 5: 0+3 nsubtype=0

Step 2&3 (start Matlab and run EpiProfile)

#### Start and Run



#### Results



# **D.** Trouble shooting

- 1. The folder of 'EpiProfile' cannot be put under a path contains a space (e.g. C:/my folder), through the data path can contain a space.
- 2. If Xcalibur is installed but EpiProfile cannot use Xcalibur to convert RAW to MS1 and MS2, then MSFileReader need be installed, which can be downloaded from https://thermo.flexnetoperations.com.
- 3. In addition to MATLAB, EpiProfile2.0 uses the following 3 toolboxes: Statistics and Machine Learning Toolbox, Curve Fitting Toolbox, and Bioinformatics Toolbox.
- 4. The application file 'xtract.exe' in the folder of 'EpiProfile' will expire on the end day of each year. Therefore, it needs be updated in the beginning of each year.